

Original Article J Biomed Transl Res 2024;25(2):53-67

ј вютеа 1ransi кез 2024;25(2):53-67 https://doi.org/10.12729/jbtr.2024.25.2.53



Received: May 23, 2024 Revised: Jun 5, 2024 Accepted: Jun 7, 2024

*Corresponding author

Byung-Suk Jeon Toxicological Evaluation Laboratory, Animal and Plant Quarantine Agency, Gimcheon 39660, Korea Tel: +82-54-912-0573 E-mail: bsjeon79@korea.kr

Copyright © 2024 Research Institute of Veterinary Medicine, Chungbuk National University. This is an Open Access article distributed under the terms of the Creative Commons Attribution Non-Commercial License (http://creativecommons.org/licenses/ by-nc/4.0/) which permits unrestricted non-commercial use, distribution, and reproduction in any medium, provided the original work is properly cited.

ORCID

Nam-Ju Kim https://orcid.org/0009-0000-4503-312X Ji-Hyun Bang https://orcid.org/0000-0003-2212-7056 Hee Yi https://orcid.org/0000-0003-0684-1882 Hyun-Ok Ku https://orcid.org/0000-0002-7890-7919 Joong-Sun Kim https://orcid.org/0000-0002-7890-7919 Ji-Yeon Kim https://orcid.org/0000-0002-5474-5951 Byung-Suk Jeon https://orcid.org/0000-0002-7981-9507

Conflict of Interest

No potential conflict of interest relevant to this article was reported.

Acknowledgements

This work was supported by research funds (B-1543073-2022-24-01) from Animal and Plant Quarantine Agency, Korea.

Comparison of *in vitro* models for drug-induced liver injury assessment

Nam-Ju Kim¹, Ji-Hyun Bang¹, Hee Yi¹, Hyun-Ok Ku¹, Joong-Sun Kim², Ji-Yeon Kim¹, Byung-Suk Jeon^{1*}

¹Toxicological Evaluation Laboratory, Animal and Plant Quarantine Agency, Gimcheon 39660, Korea ²College of Veterinary Medicine and BK21 Plus Project Team, Chonnam National University, Gwangju 61186, Korea

Abstract

Drug-induced liver injury (DILI) is considered to be a significant cause of drug wastage. To mitigate clinical DILI risks, assessing drugs using human liver models is crucial since animal studies may fall short due to species-specific liver pathway variations. Cell-based preclinical hepatotoxicity testing is often pertinent. In the present study, cells from a human liver cancer line (HepG2 and HepaRG) were cultured in both formats of 2D and 3D spheroids to explore their responses to drugs. Liver-specific marker expressions across cell lines and culture formats were also examined to assess disparities in DILI marker expressions. After treating each cell with the drugs, cytotoxicity and liver injury markers aspartate aminotransferase and alanine aminotransferase were increased. In addition, liver specific markers albumin and urea decreased in a drug concentration-dependent manner. These findings were consistent with drug sensitivity. Additionally, mRNA expression levels of cytochrome P450 enzymes (CYPs) involved in hepatocellular drug metabolism were compared following treatment with enzyme inducers. CYP1A2 and CYP2C9 were not epxressed in HepG2 cells. HepaRG cells exhibited significantly increased expression of CYP1A2, 2C9, and 3A4 post-treatment. Notably, enzyme expression was notably higher in 3D cultures than in 2D cultures. Collectively, these findings suggest that HepaRG cells and 3D cultures hold promise for evaluating DILI during early-stage drug development.

Keywords: 2D and 3D cell culture; Hep G2 cell; HepaRG cell; drug-induced liver injury; cytochrome P450

INTRODUCTION

Drug-induced liver injury (DILI) is a leading cause of preclinical or clinical drug wastage and acute liver failure [1]. Nearly 100 marketed drugs can potentially induce various types of liver injury, ranging from cell necrosis to fibrosis, or a combination thereof [2]. The severity of DILI can fluctuate based on factors such as the patient's age, gender, genetics, and underlying medical conditions. Unfortunately, the demands by regulatory agencies such as the Food and Drug Administration (FDA) for live animal testing during preclinical drug development can only detect less than half of human DILIs mainly due to species-specific differences in drug pathways and an impossibility to accurately replicate human genetics and disease conditions [3]. Challenges intrinsic to an

Ethics Approval Not applicable.

imal-based drug screening have underscored the increasing importance of in vitro human liver cultures over the past decade or so, with in vitro human liver culture becoming a cornerstone of drug discovery [4]. The process of drug discovery has historically relied on both in vitro and in vivo methods to assess the safety of compounds, with recent emphasis on high-throughput and cost-effective approaches such as cell-based screening assays. Traditionally, these in vitro cell-based assays have been performed using 2D monolayer cell culture, which are a poor mimic of in vivo conditions. Many biological functions, such as cell-cell and cell-substrate contact, are lacking in 2D monolayer cultures. This changes cell signaling pathways and reduces or deficiency of enzymes involved in metabolizing xenobiotics [5, 6]. By imitating the *in vivo* organization of natural tissues and organs, it is increasingly recognized that cells culture in 3D environments more closely mimic normal cell fuction [7–9]. In addition, many studies have shown that hepatocytes cultured in 3D form have similar structure and function to liver [10]. As a result, a great deal of effort has gone into the development of a variety of 3D cell models that hold great promise for applications in drug discovery, stem cell research, safety studies and many other cell-based analyses, bridging the gap between traditional 2D monolayer cell culture models and animal models. HepG2 cells remain popular due to their ease of handling and widespread use as one of the oldest and most well-known hepatic cell lines [11]. Previous studies have shown that culturing HepG2 cells for at least 72 hours is sufficient metabolic capacity for *in vitro* toxicity studies [12, 13]. However, HepG2 cells are known to have low activities of several metabolic enzymes and low expression levels of genes encoding cytochrome P450 enzymes (CYPs) that are essential in drug metabolism [12]. On the other hand, HepaRG cells derived from human hepatoma cell express functional liver markers including various CYPs and phase II enzymes and nuclear receptors. Additionally, HepaRG cells have been shown to be effective surrogates for primary human hepatocytes (PHHs), particularly in CYP3A4 induction studies [14]. In this study, we cultured HepG2 and HepaRG cells and exposed them to hepatotoxicants amiodarone HCl (ADR) and acetaminophen (AAP) to compare their responses. Moreover, we cultured cells in both formats of 2D and 3D spheroids to compare cytotoxicity and hepatic characteristics following drug treatment, aiming to identify which cells and culture types would be the most suitable for drug screening studies.

MATERIALS AND METHODS

Chemicals

The ADR and AAP were purchased from Sigma-Aldrich (St. Louis, MO, USA). Hepatotoxicants were chosen based on published information on hepatotoxic potential and mechanism of action [15–17]. CYP inducers Omeprazole and Rifampin were purchased from Sigma-Aldrich [18]. A concentrated stock solution of dimethylsulfoxide (DMSO, Sigma-Aldrich) was used to dissolve all chemicals and diluted with cell culture medium. Medium containing a final concentration of 0.1% DMSO was used for the untreated control group.

Cell culture and cytochrome P450 enzyme induction

The human hepatocellular carcinoma cell line (HepG2) was purchased from ATCC (Manassas, VA, USA). HepG2 cells were cultured in Eagle's minimum essential medium (EMEM, ATCC) supplemented with 10% fetal bovine serum (FBS, Merk Millipore, Burlington, MA, USA) and 1% penicillin and streptomycin (Merk Millipore). HepG2 cells were cultured in T-75 cell culture flask (Corning, New York, NY, USA) at 37°C in a 5% CO₂ humidified incubator. Cryopreserved HepaRG cells (Gibco, Braunschweig, Germany) were purchased from Thermo Fisher Scientific and maintained according to the manufacturer's instructions. HepaRG cells were cultured in William's E medium (Gibco, Billings, MT, USA) added with HepaRG[™] Thaw, Plate & General Purpose Medium supplement (Gibco) and 1% penicillin and streptomycin (Merk Millipore). HepaRG cells were cultured in T-75 cell culture flasks (Corning) at 37°C in a 5% CO₂ humidified incubator. When HepG2 and HepaRG cells reached about 80% confluency, cells were trypsinized (0.25% Trypsin-EDTA, Gibco) for 5 min, added with medium to neutralize trypsin, and centrifuged at $218 \times g$ for 5 min at 4°C. The supernatant was removed and cells were re-suspended in medium. These cells were seeded at a density of 1.5×10^4 cells/well in a total volume of 100 μ L in ultra-low attachment 96-well plates (BD Bioscience, Franklin Lakes, NJ, USA) and incubated at 37°C in a 5% CO₂ humidified incubator. For CYP enzyme induction, cells were cultured in ultra-low attachment 96-well plates. Omeprazole was used as a CYP1A2 inducer. Rifampicin was used as an inducer for CYP2C9 and CYP3A4. When cell density reached more than 80%, cells were incubated according to the condition listed in Table 1 to induce the expression of CYP enzymes.

In vitro cell viability assay

CellTiter-Glo® 2D Cell Viability Assay and CellTiter-Glo® 3D Cell Viability Assay (Promega, Madison, WI, USA) were used for the measurement of ATP as an indicator of viability, producing a luminescent output which was reported to be more sensitive than other methods [19, 20]. These assays lyse cells using proprietary lytic component. Promega has shown that the 3D cell viability assay kit has a greater ability to lyse cells in spheroids than the 2D cell viability assay kit. Therefore, we hypothesized that the 3D cell viability assay kit would be more effective in our studies. In brief, in 2D assays, 50 μ L of PBS was added to each well followed by the adding of 50 μ L CellTiter-Glo reagent which was prepared according to the manufacturer's protocol. To induce cell lysis, there contents were mixed for 2 min on a plate shaker. The samples were then incubated at room temperature for 10 min. Lysate were then transferred to a white opaque 96 well plate. Luminescent was measured using a luminometer (FlexStation III; Molecular Device, San Jose, CA, USA). For 3D assays, the plates were mixed strongly for 5 min to induce cell lysis and then incubated for 25 min at room temperature

Table 1. CYP enzyme induction methods

CYP enzyme	Inducer	Final concentration (µM)	Incubation time
CYP1A2	Omeprazole	0, 25, 50, 100	
CYP2C9 CYP3A4	Rifampin	0, 10, 25, 50	48 hr

ture to allow the luminescence signal to stabilize. Finally, $50 \ \mu\text{L}$ of the solution was transferred to a white opaque 96-well plate to be measured with a luminometer as a relative luminescence unit with an integration time of 1 sec.

Lactate dehydrogenase, aspartate aminotransferase and alanine aminotransferase assays for detecting liver injury markers

To measure enzyme leaking due to cellular damage, lactate dehydrogenase (LDH) activities were estimated using a LDH cytotoxicity detection kit (Takara Bio, Otsu, Japan) according to the manufacturer's instruction. In brief, HepG2 and HepaRG cells were seeded on 96-well cell culture plates (1.5×10^4 cells/mL) with culture medium and exposed to hepatotoxicants. Following the 24 hr treatment period, it was accordance with the manufacturer's instructions for use. The absorbance was read on a luminometer (FlexStation III) at 490 nm. Cell-free wells were used to eliminate background absorption from nanoparticles. The serum activity of aspartate aminotransferase (AST) and alanine aminotransferase (ALT) were measured using commercial kits specific for each assay. The manufacturer's instructions were followed for the protocol. ALT levels were measured using ALT activity assay kit (Sigma-Aldrich). To the cell free wells containing the reaction mixture, serum samples were added. AST levels were measured using AST activity assay kit (Sigma-Aldrich). After incubated with reaction mix for 5 min at 37°C. Absorbance was read at a wavelength of 570 nm for ALT and 450 nm for AST on a luminometer (FlexStation III). Samples were measured for absorbance at 570 nm and 450 nm every 5 min for up to 30 min. After completing the measurement, the values were calculated by subtracting the starting measurement from the last measurement.

Albumin and urea secretion

HepG2 and HepaRG cells were seeded into 96-well cell culture plates $(1.5 \times 10^4 \text{ cells/mL})$ with culture medium. On the following day, cells were treated with serial concentrations of ADR and AAP for 24 hr. To remove floating cells, the supernatants collected from each well were centrifuged at 3,000 × g for 5 min. Supernatants were then stored at -20° C until assay. For the estimation of albumin (ALB) and urea, respectively, an ALB Human ELISA kit (Abnova, Walnut, CA, USA) and a urea assay kit (Cell Biolabs, San Diego, CA, USA) were performed following the instructions of the respective supplier. Absorbance was read at a wavelength of 450 nm for ALB and 630 nm for urea on a luminometer (FlexStation III). ALB and urea levels were calculated using each standard curve. They were normalized according to cell count (1.5×10^4 cell/mL).

Qauntitative real-time polymerase chain reaction analysis of cytochrome P450 enzyme mRNA

The RNeasy Mini Kit (Qiagen, Hilden, Germany) was used to extract RNA from the cells. The concentration and purity of RNAs were assessed with an ND-1000 spectrophotometer (NanoDrop Technologies, Wilmington, DE, USA). cDNA synthesis was performed using 1 µg of the total RNA with an Omniscript RT kit (Qiageny). Polymerase chain reaction (PCR) cycling were denaturation at 95°C for 10 sec, annealing at 45°C to 64°C for 10 sec, and extension at 72°C for 20 sec for up to 40 cycles. After reverse transcription of RNA into cDNA, expression levels of CYP1A2, CYP2C9, and CYP3A4 were measured by qPCR. A CFX Real-Time PCR Detection System (Bio-Rad Laboratories, Hercules, CA, USA) was used to amplify target gene using primers. The mRNA level of each test gene (CYP1A2, CYP2C9, and CYP3A4) was normalized to β -actin according to the following formula: \triangle CT = CT (test gene) – CT (β -actin). Thereafter, relative mRNA level of each gene was calculated using the \triangle \triangle CT method: \triangle \triangle CT (test gene) = \triangle CT (test gene in the treated sample) – \triangle CT (test gene in the untreated sample). The fold change of mRNA level was expressed as the relative expression of 2^{- \triangle \triangle CT}. Induction was defined as more than 2-fold over the vehicle control and 20% or more over the positive control. The sequences of primer sets are shown in Table 2.

Data analysis

GraphPad Prism software version 5.0 was used to analyze the experimental data by towway analysis of variance (ANOVA) followed by Bonferroni post-test. The data were considered to be significant at p<0.05 (*). GraphPad Prism software 5.0 was used to calculate the 50% inhibitory concentration (IC₅₀) value (logarithmic transformation of X-values and nonlinear regression sigmoidal dose-response analysis with variable slope-with bottom and top constraints set at 0 and 100, resp.). Values are presented with \pm 95% confidence intervals.

RESULTS

Cytotoxicities of amiodarone HCI and acetaminophen to HepG2 and HepaRG cells depending on culture type

Cytotoxicities of ADR and AAP were evaluated based on cell viability using two different cell culture types (2D and 3D cultures) and compared with those of HepG2 and HepaRG cells without ADR or AAP treatment. ADR and AAP decreased cell viability of both cell lines and culture types in a dose-dependent manner (Fig. 1A and B). However, different responses were seen depending on the cell culture type. For both cells, cell viability was rapidly decreased in 2D culture, showing significant differences from 3D culture. The IC₅₀ value of ADR was found to be 44.85 μ M for HepG2 2D, 65.47 μ M for HepRG 2D, 55.64 μ M for HepG2 3D,

Prin	ner sets	Sequences
CYP1A2	Sense	5'-CACCATGGCATTGTCCCAGTCTG-3'
	Antisense	5'-CCGCGGTCAGTTGATGGAGAAGCGC-3'
CYP2C9	Sense	5'-CAGATCTGCAATAATTTTTCTC-3'
	Antisense	5'-CTTTCAATAGTAAATTCAGATG-3'
CYP3A4	Sense	5'-CTCTCATCCCAGACTTGGCCA-3'
	Antisense	5'-ACAGGCTGTTGACCATCATAAAAG-3'
β-actin	Sense	5'-CAAAGGCGAGGCTCTGTG-3'
	Antisense	5'-CCGAAAGTTGCCTTTTATGG-3'

Table 2. Primer sequences



Fig. 1. Cytotoxicities of ADR and AAP to HepG2 and HepaRG cells. HepG2 and HepaRG cells treated with ADR (0, 10, 17.8, 31.7, 56.4, 100 μ M) and AAP (0, 1, 3, 10, 30 mM) for 24 hr were incubated for 2D and 3D cultures and then used to compare cytotoxicities of ADR and AAP. Cell viability was estimated via CellTiter-Glo® luminescence assay after treatment with serial concentrations of ADR (A) and AAP (B). Values are presented as mean ± S.D. of three replicates (n = 3). p<0.05; p<0.01 (significantly different from the control). ADR, amiodarone HCI; AAP, acetaminophen.

and 109.38 μ M for HepaRG 3D. The IC₅₀ value of AAP was determined to be 14.95 mM for HepG2 2D, 7.11 mM for HepaRG 2D, 48.43 mM for HepG2 3D, and 26.19 mM for HepaRG 3D (Table 3). Overall, 2D cultures were found to be more sensitive to drugs than 3D cultures.

Detection of drug-induced hepatic injury markers

To measure the expression of LDH, AST and ALT (liver injury markers), HepG2 cells and HepaRG cells were cultured in 2D and 3D and then treated with hepatotoxicants ADR and AAP. Activities of LDH (Fig. 2A and B), AST (Fig. 2C and D), and ALT (Fig. 2E and F) were increased by treatment with ADR or AAP in a dose-dependent manner in both cell lines and culture types. However, different responses were seen depending on the type of cell culture. It was confirmed that their activities in both cells were increased more in 2D culture than in 3D culture in a drug concentration-dependent manner. HepG2 cells showed higher expression of liver injury factors after ADR treatment than HepaRG cells and 2D cultures showed higher expression of these factors than 3D cultures. HepaRG cells showed higher expression of these factors after AAP treatment than HepG2 cells and 2D cultures showed higher expression of factors than 3D cultures. These results were consistent with results of the cell viability assay.

Analysis of hepatocyte specific markers albumin and urea secretion

To determine effects of cell culture type and treatment with hepatotoxicants on expression levels of liver specific markers, differences in ALB and urea secretion between 2D and 3D cul-

	Table	3.	IC50	values	by	cell	culture	type
--	-------	----	------	--------	----	------	---------	------

IC ₅₀ value	HepG2 2D	HepaRG 2D	HepG2 3D	HepaRG 3D
Amiodarone HCl IC₅₀ (µM)	44.85	65.47	55.64	109.38
Acetaminophen IC ₅₀ (mM)	14.95	7.11	48.43	26.19

IC₅₀, 50% inhibitory concentration.



Fig. 2. Confirmation of *in vitro* drug-induced hepatotoxicity at 24 hr post-treatment. Cellular damage was estimated using LDH activities after treatment with ADR (C) and AAP (D). To compare liver injury markers ALT and AST, HepG2 and HepaRG cells treated with ADR (0, 10, 17.8, 31.7, 56.4, 100 μ M; A) and AAP (0, 1, 3, 10, 30 mM; B) for 24 hr were cultured in 2D and 3D. Values are presented as mean ± S.D. of three replicates (n = 3). *p*<0.05; *p*<0.01 (significantly different from the control). LDH, lactate dehydrogenase; AST, aspartate aminotransferase; ALT, alanine aminotransferase; ADR, amiodarone HCI; AAP, acetaminophen.

tures were analyzed. Each result was based on the control of HepG2 2D culture. Quantitative analysis of ALB showed that its expression levels were higher in 3D cultures and HepaRG cells than in 2D cultures and HepG2 cells, respectively. In addition, after culturing each cell, cells were treated with ADR or AAP to analyze changes in ALB and urea secretion depending

on the drug concentration. As a result, it was found that ALB secretion was reduced depending on the concentration of each drug (Fig. 3A and B). Urea secretion analysis gave the same results as ALB secretion (Fig. 3C and D).

Compare expression levels of cytochrome P450 enzymes after inducers treatment

We treated each cell with a CYP enzyme inducer and then examined gene expression level of the enzyme to determine which cells and culture types had higher expression levels of the CYP450 enzyme known to be involved in drug metabolism in the liver. Each result was based on the control of HepG2 2D culture. Expression levels of CYP1A2, CYP2C9, and CYP3A4 enzymes in HepG2 and HepaRG cells were found to show significant differences after treatment with drugs. In HepaRG cells, expression levels of CYP1A2, CYP2C9, and CYP3A4 were increased after treatment with inducers. Compared with HepaRG cells, HepG2 cells had relatively lower levels of CYP1A2, CYP2C9, and CYP3A4. Additionally, significant differences were seen between 2D and 3D cultures. There was no difference in the expression of CYP1A2 and CYP2C9 in 2D and 3D cultures of HepG2 cells. However, CYP3A4 had higher expression levels in 3D cultures than in 2D cultures. In HepaRG cells, all three enzymes CY-





P1A2, CYP2C9, and CYP3A4 were found to have higher expression levels in 3D cultures than in 2D cultures. Activities of CYP1A2 in 2D and 3D cultures of HepaRG were significantly (p<0.05) increased after treatment with 25 µM of omeprazole (Fig. 4A). Activities of CYP2C9 and CYP3A4 were significantly (p<0.05) increased after treatment with 10 µM rifampin (Fig. 4B and C). In HepG2 cell 2D and 3D cultures, the activity of CYP3A4 was significantly (p<0.05) increased after treatment with rifampin at 10 µM (Fig. 4C). But, CYP1A2 and CYP2C9 were not increased (p<0.05) in 2D and 3D cultures of HepG2 cells (Fig. 4A and B).

DISCUSSION

Although the use of PHHs [21], their high-throughput in vitro drug screening is severely limited by their short lifespan, high cost, limited availability and observed inter donor variability. Therefore, the immortalized hepatic lines HepG2 and HepaRG were used in this study. While HepG2 cells have many limitations in 2D culture, including the loss of hepatic-specific structure and functionality [21], the ease of handling compared to PHHs, coupled with a stable phenotype and lack of donor variability [21, 22], offers many advantages for standard operating procedures development in preclinical drug screening. On the other hand, culturing HepG2 cells in a 3D environment has been shown to restore some of the characteristics lost when cultured in a 2D environment. In addition, 3D cultured HepG2 cells has a reduced proliferation rate, self-organize and differentiate into 3D spheroids that restore lost liver structure and function [23, 24]. HepaRG cells are a hepatocellular carcinoma-derived human cell line expressing functional liver markers such as several CYPs and phase II enzymes and nuclear receptors. In CYP3A4 induction studies, HepaRG cells have been shown to be an excellent surrogate for PHHs. Therefore, we cultured cells in 2D and 3D cell culture type using HepG2 and HepaRG cells, which are most commonly used for hepatotoxicity studies. We also compared the response to drug treatment by each cell and culture type. Cell viability was measured for both 2D and 3D spheroid viability using CellTiter-Glo® ATP quantification. In parallel with the as-



Fig. 4. Comparison of CYP enzyme mRNA expression in HepG2 and HepaRG cells with 2D and 3D cultures. HepG2 and HepaRG cells treated with omeprazole (0, 25, 50, 100 μ M) and rifampin (0, 10, 25, 50 μ M) for 48 hr were cultured in 2D or 3D type. Expression levels of CYP1A2 (A), CYP2C9 (B), and CYP3A4 (C) mRNAs were evaluated after treatment with omeprazole or rifampin. Concentration-dependent responses such as increases of mRNA levels for each cell types were compared to the untreated group (DMSO) of HepG2 2D control. Values are expressed as mean ± S.D. of three replicates (n = 3). p<0.05; p<0.01 (significantly different from the control). CYP, cytochrome P450 enzyme; DMSO, dimethylsulfoxide.

sessment of viability by quantification of ATP, the LDH assays were used to determine the cytotoxicity induced by ADR and AAP. Because the LDH assay is supernatant based, the cytosolic enzyme LDH is released into the culture supernatant by compromised plasma membrane, avoiding the problems of dissociating the spheroid structure required by the other membrane consistency methods such as neutral red. When evaluated using these assays, we observed decreased viability or increased cytotoxicity induction after treatment with ADR and AAP in both culture types. When comparing the responsiveness to the drugs in each cell type, HepG2 cells were more sensitive to ADR treatment, while HepaRG cells were more sensitive to AAP treatment. Changes in cell viability based on cell culture type were also observed. For both HepG2 and HepaRG cells, it was observed that 2D cultures were more sensitive to drugs than 3D cultures. These results were also observed as a result of cytotoxicity induction. The expression of LDH increased in a drug concentration-dependent manner, with 2D cultures being more sensitive than 3D cultures. AST and ALT are markers of liver injury that increase when cells are damaged [25]. Two well-known markers of liver function are ALB and urea. ALB is a major protein synthesized in the liver and plays an important role in maintaining tumor pressure, transporting various molecules and regulating pH balance. Urea is a waste product formed in the liver during protein metabolism. The liver's ability to metabolize ammonia into urea is crucial for maintaining nitrogen balance in the body [26, 27]. HepG2 and HepaRG cells were cultured in 2D and 3D forms and each cell was treated with ADR and AAP to determine the extent of liver damage. Similar to cell viability and toxicity tests, the expression of AST and ALT increased in a drug concentration-dependent manner. Similar results were observed for the secretion of ALB and urea. Secretion decreased in a concentration-dependent manner and expression was significantly reduced in 2D cultures. In addition, HepaRG cells had higher basal expression compared to HepG2 cells when tested under the same conditions. In the same cells, expression of ALB and urea was higher in cells cultured in a 3D cultures. The hepatotoxicants ADR and AAP cause mitochondrial dysfunction [28] and liver failure [29]. ADR and AAP has been reported to cause hepatotoxicity and decrease cell viability in HepG2 [30, 31] and HepaRG [32, 33] cells. These findings are consistent with our results which showed decreased cell viability, increased expression of liver injury markers, and decreased liver function after drug treatment. The evaluation of DILI using hepatocytes relies on an accurate understanding of the characteristics of the cells and drugs, which is difficult to achieve. There are studies that show different liver cells have different responses to the same drug [34]. Alternatively, studies have shown that when a single cell is treated with different drugs and shows varying sensitivity to each drug [17]. Therefore, it is important to study liver toxicity using a variety of cells and drugs, as we have done in our study. A number of critical enzymes play an important role in cellular metabolism under disease conditions or in the metabolism of drugs. CYPs have been extensively studied since their identification because of their involvement in the metabolism of a wide range of substrates, including components of biological processes and therapeutic agents, in the human body [31, 35]. In particular, it is reported that CYPs catalyze over 95% of oxidative and reductive reaction [31]. They activate or inactivate exogenous substrates such as drugs, food additives, endogenous compounds like

steroids, xenobiotics and fatty acids [36-38]. In particular, CYP1A2, CYP2C9 and CYP3A4 enzymes are known to be the most involved in drug metabolism in human hepatocytes. CY-P1A2 is a member of the CYP family primarily found in the liver, although it is also present in other tissues. CYP1A2 is responsible for the metabolism of various substrates, including caffeine, theophylline, certain drugs, and environmental toxins like polycyclic aromatic hydrocarbons found in cigarette smoke. CYP1A2 has been extensively studied due to its role in drug metabolism and its potential impact on drug efficacy and toxicity [39]. CYP2C9 is another important CYP predominantly expressed in the liver, although it is also found in other tissues. CYP2C9 substrate specificity and catalytic properties have been extensively studied to understand its role in drug metabolism and pharmacokinetics [40]. CYP3A4 is one of the most abundant CYPs in the liver and intestine, playing a pivotal role in the metabolism of a wide range of drugs. It metabolizes approximately 50% of clinically used drugs, including statins, immunosuppressants, antivirals, and many others. Due to its broad substrate specificity and significant contribution to drug metabolism, CYP3A4 is a major determinant of drug clearance and efficacy [41]. In this study, we compared the mRNA expression of CYP enzymes in HepG2 and HepaRG cells in 2D and 3D cultures. The expression of each enzyme was induced by treating the cells with known standard inducers. HepG2 cells showed no expression of CY-P1A2 and CYP2C9 enzymes in both 2D and 3D cultures. On the other hand, the expression of CYP3A4 was increased by inducer treatment. It also showed higher expression in 3D culture. As demonstrated in various studies, HepG2 cells normally express CYP1A2 [42]. However, there is evidence that while caffeine induces CYP1A2 expression in rat hepatocytes, it may not have the same effect in human HepG2 cells [43]. HepG2 cells generally express CYP1A2, but the expression level and regulatory factors can vary depending on the experimental conditions and stimuli. There is also evidence that HepG2 cells have low expression of CYP2C9, an important enzyme involved in drug metabolism [44-46]. This low expression affects their suitability as a model for studying hepatotoxicity and drug metabolism [47, 48]. It is important to consider this limitation when using HepG2 cells as a model system for drug metabolism research. In contrast to HepG2 cells, HepaRG cells showed expression of CYP1A2, CYP2C9 and CYP3A4 enzymes in both 2D and 3D cultures. Furthermore, we found that all three enzymes were more highly expressed in 3D cultures than in 2D cultures. Interestingly, the expression of all enzymes was much higher than in HepG2 cells. HepaRG cells are a widely used model for drug metabolism studies due to their ability to express a variety of CYPs [49]. They respond to typical inducers by inducing a variety of CYPs, including CYP1A2, CY-P2C9, and CYP3A4 enzymes, which are essentially involved in human drug metabolism [50]. HepaRG cells are known to be metabolically similar to human adult hepatocytes, making them useful for drug metabolism and toxicity studies [51]. These findings corroborate our results that HepaRG cells have a higher capacity to induce CYP enzymes compared to HepG2 cells. Furthermore, our results are consistent with studies showing increased mRNA levels of CYP1A2, CYP3A4 [52, 53] and increased expression of drug metabolizing enzymes including CYP2C9 in HepaRG 3D cultures [54]. In conclusion, the development of in vitro models that can predict DILI at the preclinical stage for drug development is becoming increasingly

important. Therefore, the selection of *in vitro* models that can study a clearer and more diverse range of drug responses is important. In this study, HepG2 cells and HepaRG cells, which are often used in drug screening were compared, and each cell was cultured in 2D and 3D forms to evaluate them as *in vitro* models for drug development. The results of the study showed that the sensitivity of the cells varied depending on the type of cell, culture type and type of drug. In addition, HepaRG cells with higher expression of CYP enzymes involved in human drug metabolism are more suitable for drug metabolism evaluation, and 3D culture is more suitable as a model for drug metabolism evaluation than 2D culture. However, the results of CYP enzymes need to be validated by further experiments on metabolic mechanisms through enzyme activity measurement or liquid chromatography-tandem mass spectrometry metabolite.

REFERENCES

- 1. Kaplowitz N. Idiosyncratic drug hepatotoxicity. Nat Rev Drug Discov 2005;4:489-499.
- 2. Abboud G, Kaplowitz N. Drug-induced liver injury. Drug Saf 2007;30:277-294.
- Olson H, Betton G, Robinson D, Thomas K, Monro A, Kolaja G, Lilly P, Sanders J, Sipes G, Bracken W, Dorato M, Van Deun K, Smith P, Berger B, Heller A. Concordance of the toxicity of pharmaceuticals in humans and in animals. Regul Toxicol Pharmacol 2000;32:56-67.
- Lin C, Khetani SR. Advances in engineered liver models for investigating drug-induced liver injury. BioMed Res Int 2016;2016:1829148.
- Edmondson R, Broglie JJ, Adcock AF, Yang L. Three-dimensional cell culture systems and their applications in drug discovery and cell-based biosensors. Assay Drug Dev Technol 2014;12:207-218.
- Wrzesinski K, Fey SJ. From 2D to 3D: a new dimension for modelling the effect of natural products on human tissue. Curr Pharm Des 2015;21:5605-5616.
- Zhang X, Yang ST. High-throughput 3-D cell-based proliferation and cytotoxicity assays for drug screening and bioprocess development. J Biotechnol 2011;151:186-193.
- Fey SJ, Wrzesinski K. Determination of drug toxicity using 3D spheroids constructed from an immortal human hepatocyte cell line. Toxicol Sci 2012;127:403-411.
- Gunness P, Mueller D, Shevchenko V, Heinzle E, Ingelman-Sundberg M, Noor F. 3D organotypic cultures of human HepaRG cells: a tool for *in vitro* toxicity studies. Toxicol Sci 2013;133:67-78.
- 10. Bell CC, Hendriks DFG, Moro SML, Ellis E, Walsh J, Renblom A, Fredriksson Puigvert L, Dankers ACA, Jacobs F, Snoeys J, Sison-Young RL, Jenkins RE, Nordling Å, Mkrtchian S, Kevin Park B, Kitteringham NR, Goldring CEP, Lauschke VM, Ingelman-Sundberg M. Characterization of primary human hepatocyte spheroids as a model system for drug-induced liver injury, liver function and disease. Sci Rep 2016;6:25187.
- Hart SN, Li Y, Nakamoto K, Subileau E, Steen D, Zhong X. A comparison of whole genome gene expression profiles of HepaRG cells and HepG2 cells to primary human hepatocytes and human liver tissues. Drug Metab Dispos 2010;38:988-994.
- 12. O'Brien PJ, Edvardsson A. Validation of a multiparametric, high-content-screening assay for

predictive/investigative cytotoxicity: evidence from technology transfer studies and literature review. Chem Res Toxicol 2017;30:804-829.

- O'Brien PJ, Irwin W, Diaz D, Howard-Cofield E, Krejsa CM, Slaughter MR, Gao B, Kaludercic N, Angeline A, Bernardi P, Brain P, Hougham C. High concordance of drug-induced human hepatotoxicity with *in vitro* cytotoxicity measured in a novel cell-based model using high content screening. Arch Toxicol 2006;80:580-604.
- McGinnity DF, Zhang G, Kenny JR, Hamilton GA, Otmani S, Stams KR, Haney S, Brassil P, Stresser DM, Riley RJ. Evaluation of multiple *in vitro* systems for assessment of CYP3A4 induction in drug discovery: human hepatocytes, pregnane X receptor reporter gene, and Fa2N-4 and HepaRG cells. Drug Metab Dispos 2009;37:1259-1268.
- Gomez-Lechon MJ, Lahoz A, Gombau L, Castell JV, Donato MT. *In vitro* evaluation of potential hepatotoxicity induced by drugs. Curr Pharma Des 2010;16:1963-1977.
- Tolosa L, Pinto S, Teresa Donato M, Lahoz A, Castell JV, Enrique O'Connor J, José Gómez-Lechón M. Development of a multiparametric cell-based protocol to screen and classify the hepatotoxicity potential of drugs. Toxicol Sci 2012;127:187-198.
- Tham NTT, Hwang SR, Bang JH, Yi H, Park YI, Kang SJ, Kang HG, Kim YS, Ku HO. High-content analysis of *in vitro* hepatocyte injury induced by various hepatotoxicants. J Vet Sci 2019;20:34-42.
- Huang SM, Temple R, Throckmorton DC, Lesko LJ. Drug interaction studies: study design, data analysis, and implications for dosing and labeling. Clin Pharmacol Ther 2007;81:298-304.
- Kleijn A, Kloezeman JJ, Balvers RK, Kaaij M, Dirven CMF, Leenstra S, Lamfers MLM. A systematic comparison identifies an ATP-based viability assay as most suitable read-out for drug screening in glioma stem-like cells. Stem Cells Int 2016;2016:5623235.
- Crouch SPM, Kozlowski R, Slater KJ, Fletcher J. The use of ATP bioluminescence as a measure of cell proliferation and cytotoxicity. J Immunol Methods 1993;160:81-88.
- LeCluyse EL, Witek RP, Andersen ME, Powers MJ. Organotypic liver culture models: meeting current challenges in toxicity testing. Crit Rev Toxicol 2012;42:501-548.
- 22. Gomez-Lechon MJ, Donato MT, Lahoz A, Castell JV. Cell lines: a tool for *in vitro* drug metabolism studies. Curr Drug Metab 2008;9:1-11.
- Wrzesinski K, Magnone MC, Hansen LV, Kruse ME, Bergauer T, Bobadilla M, Gubler M, Mizrahi J, Zhang K, Andreasen CM. HepG2/C3A 3D spheroids exhibit stable physiological functionality for at least 24 days after recovering from trypsinisation. Toxicol Res 2013;2:163-172.
- Bokhari M, Carnachan RJ, Cameron NR, Przyborski SA. Culture of HepG2 liver cells on three dimensional polystyrene scaffolds enhances cell structure and function during toxicological challenge. J Anat 2007;211:567-576.
- McGill MR. The past and present of serum aminotransferases and the future of liver injury biomarkers. EXCLI J 2016;15:817-828.
- 26. Spinella R, Sawhney R, Jalan R. Albumin in chronic liver disease: structure, functions and therapeutic implications. Hepatol Int 2016;10:124-132.
- 27. Nagamani SCS, Ali S, Izem R, Schady D, Masand P, Shneider BL, Leung DH, Burrage LC.

Biomarkers for liver disease in urea cycle disorders. Mol Genet Metab 2021;133:148-156.

- Yan M, Huo Y, Yin S, Hu H. Mechanisms of acetaminophen-induced liver injury and its implications for therapeutic interventions. Redox Biol 2018;17:274-283.
- 29. Essrani R, Mehershahi S, Essrani RK, Ravi SJK, Bhura S, Sudhakaran A, Hossain M, Mehmood A. Amiodarone-induced acute liver injury. Case Rep Gastroenterol 2020;14:87-90.
- Behrends V, Giskeødegård GF, Bravo-Santano N, Letek M, Keun HC. Acetaminophen cytotoxicity in HepG2 cells is associated with a decoupling of glycolysis from the TCA cycle, loss of NADPH production, and suppression of anabolism. Arch Toxicol 2019;93:341-353.
- Zahno A, Brecht K, Morand R, Maseneni S, Török M, Lindinger PW, Krähenbühl S. The role of CYP3A4 in amiodarone-associated toxicity on HepG2 cells. Biochem Pharmacol 2011; 81:432-441.
- 32. Allard J, Bucher S, Massart J, Ferron PJ, Le Guillou D, Loyant R, Daniel Y, Launay Y, Buron N, Begriche K, Borgne-Sanchez A, Fromenty B. Drug-induced hepatic steatosis in absence of severe mitochondrial dysfunction in HepaRG cells: proof of multiple mechanism-based toxicity. Cell Biol Toxicol 2021;37:151-175.
- McGill MR, Yan HM, Ramachandran A, Murray GJ, Rollins DE, Jaeschke H. HepaRG cells: a human model to study mechanisms of acetaminophen hepatotoxicity. Hepatology 2010;53:974-982.
- Kang SJ, Lee HM, Park YI, Yi H, Lee H, So B, Song JY, Kang HG. Chemically induced hepatotoxicity in human stem cell-induced hepatocytes compared with primary hepatocytes and HepG2. Cell Biol Toxicol 2016;32:403-417.
- El-Sherbeni AA, El-Kadi AOS. Microsomal cytochrome P450 as a target for drug discovery and repurposing. Drug Metab Rev 2017;49:1-17.
- Capdevila JH, Harris RC, Falck JR. Microsomal cytochrome P450 and eicosanoid metabolism. Cell Mol Life Sci 2002;59:780-789.
- Nebert DW, Dalton TP. The role of cytochrome P450 enzymes in endogenous signalling pathways and environmental carcinogenesis. Nat Rev Cancer 2006;6:947-960.
- Guengerich FP. Human cytochrome P450 enzymes. In: de Montellano PRO (ed.). Cytochrome P450: structure, mechanism, and biochemistry. Cham: Springer; 2015. p. 523-785.
- Zhou SF, Wang B, Yang LP, Liu JP. Structure, function, regulation and polymorphism and the clinical significance of human cytochrome P450 1A2. Drug Metab Rev 2010;42:268-354.
- Wang B, Wang J, Huang SQ, Su HH, Zhou SF. Genetic polymorphism of the human cytochrome P450 2C9 gene and its clinical significance. Curr Drug Metab 2009;10:781-834.
- 41. Peter Guengerich F. Cytochrome P-450 3A4: regulation and role in drug metabolism. Annu Rev Pharmacol Toxicol 1999;39:1-17.
- Xie C, Pogribna M, Word B, Lyn-Cook L Jr, Lyn-Cook BD, Hammons GJ. *In vitro* analysis of factors influencing CYP1A2 expression as potential determinants of interindividual variation. Pharmacol Res Perspect 2017;5:e00299.
- Washio I, Maeda M, Sugiura C, Shiga R, Yoshida M, Nonen S, Fujio Y, Azuma J. Cigarette smoke extract induces CYP2B6 through constitutive androstane receptor in hepatocytes. Drug Metab Dispos 2011;39:1-3.

- 44. Horn JR, Hansten PD. Get to know an enzyme: CYP1A2. Pharm Times 2007;73:76.
- 45. Tolosa L, Gómez-Lechón MJ, Pérez-Cataldo G, Castell JV, Donato MT. HepG2 cells simultaneously expressing five P450 enzymes for the screening of hepatotoxicity: identification of bioactivable drugs and the potential mechanism of toxicity involved. Arch Toxicol 2013;87:1115-1127.
- Jover R, Bort R, Gómez-Lechón MJ, Castell JV. Re-expression of C/EBPα induces CYP2B6, CYP2C9 and CYP2D6 genes in HepG2 cells. FEBS Lett 1998;431:227-230.
- 47. Sekretarska J, Szczepaniak J, Sosnowska M, Grodzik M, Kutwin M, Wierzbicki M, Jaworski S, Bałaban J, Daniluk K, Sawosz E. Influence of selected carbon nanostructures on the CYP2C9 enzyme of the P450 cytochrome. Materials 2019;12:4149.
- Negoro R, Tasaka M, Deguchi S, Takayama K, Fujita T. Generation of HepG2 cells with high expression of multiple drug-metabolizing enzymes for drug discovery research using a PITCh system. Cells 2022;11:1677.
- 49. Gerets HHJ, Tilmant K, Gerin B, Chanteux H, Depelchin BO, Dhalluin S, Atienzar FA. Characterization of primary human hepatocytes, HepG2 cells, and HepaRG cells at the mRNA level and CYP activity in response to inducers and their predictivity for the detection of human hepatotoxins. Cell Biol Toxicol 2012;28:69-87.
- Kanebratt KP, Andersson TB. HepaRG cells as an *in vitro* model for evaluation of cytochrome P450 induction in humans. Drug Metab Dispos 2008;36:137-145.
- Ueyama T, Tsuji S, Sugiyama T, Tada M. Fluorometric evaluation of CYP3A4 expression using improved transgenic HepaRG cells carrying a dual-colour reporter for CYP3A4 and CY-P3A7. Sci Rep 2017;7:2874.
- Mandon M, Huet S, Dubreil E, Fessard V, Le Hégarat L. Three-dimensional HepaRG spheroids as a liver model to study human genotoxicity *in vitro* with the single cell gel electrophoresis assay. Sci Rep 2019;9:10548.
- Higuchi Y, Kawai K, Kanaki T, Yamazaki H, Chesné C, Guguen-Guillouzo C, Suemizu H. Functional polymer-dependent 3D culture accelerates the differentiation of HepaRG cells into mature hepatocytes. Hepatol Res 2016;46:1045-1057.
- Ott LM, Ramachandran K, Stehno-Bittel L. An automated multiplexed hepatotoxicity and CYP induction assay using HepaRG cells in 2D and 3D. SLAS Discov 2017;22:614-625.